

## Android / iOS Screens

This table lists each screen that needs to be developed along with a brief description of their purpose.

Copy	High-Fidelity Mockup
------	----------------------

## App Loading

- Shown when user starts app and app is loading
- The top navbar is optional
- Main panel will show a static image or an animation.

\*\*\*\*\* Can we make the Sana logo larger?\*\*\*\*\*

The image shows a dark blue background with a pattern of wavy, horizontal lines in a slightly lighter shade of blue, creating a sense of movement. In the upper center, the word "sana" is written in a white, lowercase, sans-serif font. Below it, the tagline "Relief made possible." is written in a smaller, white, sans-serif font.

sana  
Relief made possible.

# sana




Looks like you haven't connected your Sana Device to the app. Click the button below to connect your Sana now.

**Connect your Sana Device**

## Login

- Shown when the app is in a loaded state, but the user session state is inactive.
- The user needs to provide their email and password in order to login to the app.
  - If a user enters invalid account credentials then we need to use form validation to present the error to the user.
- Alternatively they can create an account which takes them to the create account page.
- If the user needs to reset their password they can click the Password Reset button which opens a modal on their screen asking them if they want to proceed to a web page.



The image shows a mobile app login screen for 'sana'. At the top, the 'sana' logo is in a light blue font. Below it, there are two tabs: 'Log In' (highlighted with a blue underline) and 'Signup'. The main area contains two input fields: 'Email' with the placeholder 'johnsmith@example.com' and 'Password' with the placeholder 'Password'. A blue 'Show' link is next to the password field. At the bottom, there is a large blue 'Log In' button and a 'Forgot Password?' link.

sana

Log In Signup

Email

johnsmith@example.com

Password

Password Show

Log In

Forgot Password?


## Account Creation

- Shown when the user needs to create an account.
- Back button takes a user back to login screen.
- If the user clicks the Create Account button and their form is *invalid* then we need to show the user that there is an issue.
- If the user clicks the Create Account button and their form is *valid* then we create the user account and log the user into the app.

Change DOB to Birthdate

Do we need to add a 'Show' on the Confirm Password lined as well?

On Warning, You must be... Sana Device We want to ensure when we say Sana Device, the D is always capital.



The image shows a mobile app interface for the 'sana' app. At the top, the 'sana' logo is in a light blue font. Below it, there are two tabs: 'Log In' and 'Signup'. The 'Signup' tab is selected, indicated by a light blue underline. The main area contains four input fields: 'Email' with the placeholder 'johnsmith@example.com', 'Birthdate' with the placeholder '01/02/1997', 'Password' with the placeholder 'Password' and a 'Show' link to its right, and 'Confirm Password' with the placeholder 'Confirm Password'. At the bottom, there is a large blue button with the text 'Sign Up' in white.

sana

Log In Signup

Email

johnsmith@example.com

Birthdate

01/02/1997

Password

Password Show

Confirm Password

Confirm Password

Sign Up

sana

Log In

Signup

Email

johnsmith@example.com

Email Required

Birthdate

01/02/1997

Birthdate Required

Password

Password

Show

Password Required

Confirm Password

Confirm Password

Confirm Password Required

Sign Up

sana

Log In

Signup

Email

johnsmith@example.com

Email Required

Birthdate

01/02/1997

Birthdate Required

Password

Password

Show

Password Required

Confirm Password

Confirm Password

Confirm Password Required

Sign Up

sana

Log In

Signup

Email

johnsmith@example.com

Email Required

Birthdate

01/02/1997

Birthdate Required

Password

Password

Show

Password Required

Confirm Password

Confirm Password

Confirm Password Required

Sign Up

--	--

## Email Validation

- Shown after login prior to their email being verified.
- They enter a verification code which they receive via email and click submit. If the verification code is valid then they move onto the next stage of their account setup, which is giving us their profile information.
- If their verification code isn't working we give them the option of requesting a new verification code.

We sent you an email containing a validation code. Please enter the code in the field below and click Next to validate your email.

Sana Email Validation Code



## Email Validation

We sent you an email containing a validation code. Please enter the code in the field below and click Next to validate your email.

Validation Code


A9P9123

Next

[Send another Validation Code](#)

## Forgot Password Request Validation Code

- Back button takes user back to login screen
- Clicking the button takes them to the Forgot password enter validation code screen.

Forgot Password

Send me an email  
to reset my password


Email

Next

## Forgot Password Enter Validation Code

- Back button takes user back to login screen
- Reset password button attempts to reset their password via the server.

Instead of Token, please use "Validation Code"

 **Forgot Password**

We sent you an email containing a validation code. Please use the form below to reset your password.

Validation Code

A9P9123

Password

Password [Show](#)

Confirm Password

Confirm Password

Reset Password

[Send another Validation Code](#)

## 1st Logon Introduction

Welcome to your Sana!

Please take a few moments to help us get to know you.

The Sana logo is displayed in a light blue, lowercase, sans-serif font against a dark navy blue background.

Welcome to your Sana!

Please take a few moments to help us  
get to know you.

Got it

## 1st Logon Setup Profile

- Shown on first ever login to the app.
- Designed to gather profile and demographic data.
- Nothing is required, but we don't explicitly say skip as we want them to fill in the data.

Gender options:

Female

Male

Do not identify

Decline to Answer

Update DOB to Birthdate



## Create Profile

### Step 1 of 3

First Name

John

Last Name

Smith

Gender

☐

Female

☐

Male

☐

Do not identify

☐

Decline to answer

Race

☐

America Indian or Alaska Native

☐

Asian

☐

Black or African American

☐

Native Hawaiian or other Pacific islander

☐

White



## Create Profile

### Step 1 of 3

First Name

John

Last Name

Smith

Gender

☐

Female

☐

Male

☐

Do not identify

☐

Decline to answer

Race

☐

America Indian or Alaska Native

☐

Asian

☐

Black or African American

☐

Native Hawaiian or other Pacific islander

## 1st Logon Setup Data Collection

- Shown on first ever login to the app.
- Designed to gather data from user on average sleep, pain, mood, etc.
- Nothing is required, but we don't explicitly say skip as we want them to fill in the data.

(Intro) Over the past month, how would you describe:

Sleep Scale: Terrible to Excellent

Mood Scale: Depressed to Happy

 **Create Profile**

Step 2 of 3

Over the past month,  
how would you describe:

**Average Sleep Quality**  
  
Terrible Excellent

**Average Pain**  
  
Debilitating Faint

**Average Mood**  
  
Depressed Happy

Next




## Create Profile


### Step 2 of 3

Over the past month,  
how would you describe:


#### Average Sleep Quality



#### Average Pain




#### Average Mood



Next


## 1st Logon Setup Video

Click below to learn more about your Sana Device.



## Forgot Password Request Validation Code

- Back button takes user back to login screen
- Clicking the button takes them to the Forgot password enter validation code screen.

Forgot Password

Send me an email  
to reset my password

Email

Next


## Forgot Password Enter Validation Code

- Back button takes user back to login screen
- Reset password button attempts to reset their password via the server.

Obtain the Validation Code from your email, then use the form below to reset your password.

Validation Code (instead of Token)

Send another Validation Code (Instead of Send another token)

 **Forgot Password**

We sent you an email containing a validation code. Please use the form below to reset your password.

Validation Code

A9P9123

Password

Password

Show

Confirm Password

Confirm Password

Reset Password

Send another Validation Code

## Main App Layout

- App has the following layout:
  - Top Nav
  - Inner Panel, scrollable
  - Bottom Nav
- When the user logs in the My Sana tab (lower panel, left side) is active.
- When the Profile tab is clicked, the inner panel refreshes to show the Profile screen.
- When the More icon is clicked, the inner panel updates to show a menu of options, including videos, FAQ, support, feedback, guides, etc.
- Long term we're going to want to add other panels like notifications, surveys, my journal, etc.


## My Sana - Device Not Registered, with historical data

- If the user has not registered a mask then they see the following screen when the My Sana tab is active.
- The app automatically starts looking for Sana Masks broadcasting via BLE.

Welcome back 'Name' This is at the top before any verbiage

Looks like you haven't connected your Sana Device to the app. Click the button below to connect your Sana now.

(in Button) Connect your Sana Device



sana



## We cannot find your Sana Device

Please ensure your Sana Device is turned on. If the problem persists, please turn your Sana Device off and back on. If you continue to have difficulty, please contact [customer support](#).




### Last Check-In

5 minutes ago

Sleep



Pain



Mood



This Week


Week 10

--	--

## My Sana - Device Not Registered, without historical data

Looks like you haven't connected your Sana Device to the app. Click the button below to connect your Sana now.

(in Button) Connect your Sana Device



## Registration Setup Step 1

### Connection Setup (Not Registration Setup)

You are just a few moments away from accessing your relief.

The first step to connecting your Sana Device is to make sure the power switch is in the “On” position.

If you need helping locating the “on” position, check out our FAQs ([hyperlink to FAQs](#))




## Connection Setup

**You are just a few moments away  
from accessing your relief**

The first step to connecting your Sana Device is to make sure the power switch is in the "On" position.

If you need helping locating the "on" position, check out our [FAQs](#)



Next



## Registration Setup Step 2

Connection Setup (Not Registration Setup)

The Sana Device requires your phone bluetooth to be turned on.



### Connection Setup

The Sana Device requires your  
phone bluetooth to be turned on



Next



## Registration Setup Step 3

### Connection Setup (Not Registration Setup)

To connect your Sana Device with your phone you'll need the Sana Device ID which can be located on the Packaging Sleeve or Warning Card.

Below you will find all Sana Devices that are available to connect with your phone via bluetooth.

Select the option below that matches your Sana Device ID. If you don't see your Sana Device ID, [click here](#) for help.



## Connection Setup

**To connect your Sana Device with your phone you'll need the Sana Device ID which can be located on the Packaging Sleeve or Warning Card.**

Below you will find all Sana Devices that are available to connect with your phone via bluetooth.

Select the option below that matches your Sana Device ID. If you don't see your Sana Device ID, [click here](#) for help.

Sana-APL9IA

Sana-MK191K

**Sana-OIJ89A**

Sana-APL9IA

Sana-MK191K

Next




## Registration Step 3 - Connection Error

### Connection Setup (Not Registration Setup)

Sorry, there has been an issue connecting with your Sana Device. Please turn your Sana device off and then back on and try connecting again.

If you continue to have difficulty connecting your Sana Device, please contact [support](#). (link to support page)



## Registration Step 3 - No devices found

- Any masks are presented in a clickable list.
- If the user sees their mask ID then they click on it and immediately the screen updates to show the default My Sana page.

### Connection Setup (Not Registration Setup)

To connect your Sana Device with your phone you'll need the Sana Device ID which can be located on the Packaging Sleeve or Warning Card.

Below you will find all Sana Devices that are available to connect with your phone via bluetooth.

Select the option below that matches your Sana Device ID. If you don't see your Sana Device ID, [click here](#) for help.



## Connection Setup

**To connect your Sana Device with your phone you'll need the Sana Device ID which can be located on the Packaging Sleeve or Warning Card.**

Below you will find all Sana Devices that are available to connect with your phone via bluetooth.

Select the option below that matches your Sana Device ID. If you don't see your Sana Device ID, [click here](#) for help.



Searching....

## Mask Not Found Help Required


### Sana Connection Help

If you are having trouble connecting your Sana Device, please try the following.

Turn your Sana Device off and back on again.

Move your Sana Device closer to your phone.

If your problem persists, please contact customer [support](#). (link to support page)




## Registration Step 4 - All Good

### Connection Setup (Not Registration Setup)

You're all set! Thanks for connecting your Sana Device.

Click the button below to Start your Relief

Otherwise click the back button in the upper left navigation panel to return to the My Relief homepage.



## Registration step 4- Issue with account

### Connection Setup (Not Registration Setup)

There is an issue with your account that requires you to contact customer service.

Please call them at: 720-310-3370

(Button) - Restart Connection

(Button) - Cancel Connection



### Connection Setup

**There is an issue with your account that requires you to contact customer service.**

Please call them at:

720-310-3370

**Restart Connection**

**Cancel Connection**



## Connection Setup

**We are sorry, but there appears to be an issue with your account. Please contact customer service so we can get you back to Your Relief as soon as possible.**

Call us at:


720-310-3370

**Restart Connection**

**Cancel Connection**


## My Sana - Device Registered and Connected w/ No Usage


- If a mask is registered and the mask is connected via BLE then the user sees this screen.
- If you click the start a new session button you will move to the Session Prior to Start screen.




## My Sana - Device Registered and Connected

- If a mask is registered then the user sees this screen.
- If you click the start a new session button you will move to the Session Prior to Start screen.
- The history section shows a card for each prior usage. On the card we'll show the date, duration, and any pain vas scores we can gather.







## Account Issue

- Thrown when a user is trying to use a mask that has an account issue.

There is an issue with your account that requires you to contact customer service.

Please call them at: 720-310-3370

(Button) - Restart Connection



## My Relief

**There is an issue with your account that requires you to contact customer service.**

Please call them at:

720-310-3370

**Restart Connection**



## My Relief

**We are sorry, but there appears to be an issue with your account. Please contact customer service so we can get you back to Your Relief as soon as possible.**

Call us at:


**720-310-3370**

**Restart Connection**

## BLE Not On

- If bluetooth on the device (phone, tablet) is turned off then the app presents this message telling the user to turn their BLE on and return to the app.
- Once they do this message will disappear.

Your Sana Device requires your phone bluetooth to be turned on.



## My Sana - Device Registered and Not Connected

- If a device is not connected then the My Sana page shows an error message and does not let the user start a session or see the connected icons in the top navbar.
- The app will continue to look for the mask every 2 seconds until it finds it.
- When it does the screen updates to the My Sana - Device Registered and Connected screen.

### We cannot find your Sana Device


Please ensure your Sana Device is turned on. If the problem persists, please turn your Sana Device off and back on. If you continue to have difficulty, please contact [customer support](#). (hyperlink to support)



## My Relief

### We cannot find your Sana Device


Please ensure your Sana Device is turned on. If the problem persists, please turn your Sana Device off and back on. If you continue to have difficulty, please contact [customer support](#).





## My Sana - Device Battery too Low

- If a Sana mask does not have enough battery power then the user sees an error and cannot start a new session.

Your Sana Device does not have enough charge. Please charge your Sana for at least one hour and then try again.







## My Relief Optional Data Input

- The user has to either press skip or somewhere on the slider for each question.
- Continue button not active until each section has been answered.

Sleep Scale: Terrible to Excellent

Mood Scale: Depressed to Happy


The screenshot shows the 'sana' app interface with a dark background. At the top, the logo 'sana' is displayed in a light blue font. Below the logo, there are three sections, each containing a slider and a 'Skip' button:

- Average Sleep:** The slider ranges from 'Terrible' to 'Excellent'. A blue dot is positioned at the far right end, indicating 'Excellent'.
- Average Pain:** The slider ranges from 'Debilitating' to 'Faint'. A blue dot is positioned at the far left end, indicating 'Debilitating'.
- Average Mood:** The slider ranges from 'Depressed' to 'Happy'. A blue dot is positioned in the middle, indicating a neutral mood.

At the bottom of the screen, there is a large blue button labeled 'Continue'.


sana

Average Sleep



Skip

Average Pain



Skip

Average Mood




Skip

Continue

## Session Prior to Start

- The top section shows the counter and the state (one of countdown or in-progress).
- The Elapsed time counts up once the session starts. The Remaining counter counts down.
- The ? opens a modal with a brief description of how to use the interface.
- The Start button lets you start the 5 second countdown. This takes you to the Session Countdown Active screen. This is shown on the Session Stopped Early screen.
- The X in the upper left opens a modal asking if you want to save or discard the session data or cancel the modal.


Header: **My Relief Countdown**



## Session Countdown Active

- The top section shows the counter and the state (either in the countdown or the in-progress).
- The Elapsed time counts up once the session starts. The Remaining counter counts down. They haven't changed yet as we're still in the countdown phase.
- The ? opens a modal with a brief description of how to use the interface.
- The Pause button lets you pause the 5 second countdown. This takes you to the Session Countdown Paused screen. This is shown on the Session Stopped Early screen.
- The X in the upper left opens a modal asking if you want to save or discard the session data or cancel the modal.


Header: **My Relief Countdown**



## Session Countdown Paused

- The top section shows the counter and the state (either in the countdown or the in-progress).
- The Elapsed time counts up once the session starts. The Remaining counter counts down. They haven't changed yet as we're still in the countdown phase.
- The ? opens a modal with a brief description of how to use the interface.
- The Resume button lets you resume the 5 second countdown. This takes you back to the Session Countdown Active screen.
- The X in the upper left opens a modal asking if you want to save or discard the session data or cancel the modal.
- After the 5 second counter is complete the screen changes to the Session in Progress screen.


Header: **My Relief Countdown**



## Session In Progress

- The top section shows the counter and the state which is In Progress.
- The Elapsed time counts up once the session starts. The Remaining counter counts down.
- The ? opens a modal with a brief description of how to use the interface.
- The stop button or X in the upper left open a modal asking if you want to save or discard the session data or cancel the modal.
- After the 15 minute screen the Session Complete screen is shown.

Header: **My Relief In Progress**




## Session In Progress Help Screen

- 

If you are having issues with your Sana Device,  
please contact customer support via phone at:

720-310-3370


Or via email at:  
[support@sana.io](mailto:support@sana.io)



## Session Stopped Early

- If you click discard the session data will not be saved.
- If you click save the session data will be saved
- If you click cancel the modal closes.


Let's make the following words lowercase (you, like and session)



## Session Complete

- If you click discard the session data will not be saved.
- If you click save the session data will be saved
- If the user doesn't do anything for 60 seconds, then we save the session and close the modal.


Let's make the following words lowercase (you, like and session)



## Me - View Insights

Change Mood Improvement, Sleep Improvement, and Pain Improvement to Mood, Sleep and Pain.

Also, they want to see another line chart that has an x- and y-axis. No labels on either axes. Do your best, but don't spend too much time on it.



sana

Insights

Profile

Your Averages

5 minutes ago

Sleep




Pain



Mood



This week



Sleep

Pain

Mood

--	--

sana

Insights

Profile

### Your Averages


Sleep



Pain



Mood



This week



Once you have more data you will see changes in your mood, sleep and pain levels here.



My Relief



Me



More

--	--

sana

Insights

Profile

### March Usage

Su	Mo	Tu	We	Th	Fr	Sa
01	02	03	04	05	06	07
08	09	10	11	12	13	14
15	16	17	18	19	20	21
22	23	24	25	26	27	28
29	30	31				

Once you have more data you will see changes in your mood, sleep and pain levels here.



My Relief



Me



More

--	--

## Me - View Profile

Update DOB to Birthdate

sana

Insights

Profile

First Name

John

Last Name

Smith

Gender

Male

Race

Caucasian

Email

johnsmith@mail.com

Mobile Number

+1 2031045321

Address

162 Quigley Islands Suite 220

Birthdate

08/14/1978

## Me - Edit Profile

Gender options:

Female

Male

Do not identify

Decline to Answer

Update DOB to Birthdate



## Edit Profile

First Name

John

Last Name

Smith

Gender



Female



Male



Do not identify



Decline to answer

Race



America Indian or Alaska Native



Asian



Black or African American



Native Hawaiian or other Pacific islander




White

## More

- Initially each of these opens a link to a website that displays the content.
- Eventually we might want to open new tabs.
- The Logout link throws a confirm logout modal.
- If they click the Back button in the top nav then they go back to the My Sana page.
- If they click Sana Mask Details they go to the Sana Mask Details screen.

My Sana Device Details (Instead of Sana Mask Details)



## More

Product Support



FAQs



My Sana Device Details



Video



Account



Reset Password



Terms and Conditions



Logout





## Terms and Conditions

### Name

Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam, quis nostrud exercitation ullamco laboris nisi ut aliquip ex ea commodo consequat. Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint occaecat cupidatat non proident, sunt in culpa qui officia deserunt mollit anim id est laborum.

Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam, quis nostrud exercitation ullamco laboris nisi ut aliquip ex ea commodo consequat. Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint occaecat cupidatat non proident, sunt in culpa qui officia deserunt mollit anim id est laborum.


Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam, quis nostrud exercitation ullamco laboris nisi ut aliquip ex ea commodo consequat. Duis aute irure dolor in reprehenderit in voluptate velit esse cillum dolore eu fugiat nulla pariatur. Excepteur sint

## Sana Mask Details

- Just shows various details on the mask. Like hardware and firmware versions .
- Clicking back takes you back to the More screen.

My Sana Device (instead of Sana Health Mask)

Change “Unregistered Mask” to “Disconnect Sana Device”



## Unregister Mask

My Sana Device (instead of Sana Health Mask)

Click the button below to disconnect your Sana Device

(Button) : Disconnect My Sana



My Sana Device

Click the button below to  
disconnect your Sana Device

Disconnect My Sana


## Unregister Mask Confirmation

- Thrown when a user clicks unregister mask on the More screen.

My Sana Device (instead of Sana Health Mask)

Are you sure you want to disconnect?

Please confirm your decision to disconnect your Sana Device.




## Logout Confirmation

- Thrown when a user clicks the logout button on the More screen.

Are you sure you want to logout?

Please confirm your logout below



## Product Support

To contact customer support please call

720-310-3370

Or send email to:  
[support@sana.io](mailto:support@sana.io)



## Product Support

To contact customer support  
please call

720-310-3370

Or send email to:

[support@sana.io](mailto:support@sana.io)

## FAQS

To see our Frequently Asked Questions (FAQs) on our website, please click below.



## FAQ

To see our Frequently Asked Questions (FAQs) on our website, please click below.

[Browse FAQs](#)



## FAQs

To see our Frequently Asked Questions (FAQs) on our website, please click below.

[Browse FAQs](#)



## Video

Click below to learn more about your Sana Device.



## Video

Click below to learn more about  
your Sana Device.





## Video


To see helpful video content on  
our website please click below.

[Browse Videos](#)

More Account	
--------------	--

## Reset Password By Authenticated User

Button should say Reset My password

 **Reset Password**

Current Password

Current Password [Show](#)

Password

Password [Show](#)

Confirm Password

Confirm Password

[Reset My password](#)

